

PyroSmart NextGen™ 重組蛋白試劑的確效方法

Associates of Cape Cod (ACC) 為了保育鱈魚而研發的重組蛋白試劑 PyroSmart NextGen™ Recombinant Cascade Reagent (PSNG rCR) 較生物性原料具有批次間一致性佳、來源穩定、能避免葡聚糖 (1,3-β-D-Glucans) 交叉污染等等優勢，目前為檢測內毒素的替代方法學，因此終端使用者必須確認此替代方法效能 (Performance) 等同或優於傳統鱈魚血試劑 Limulus Amoebocyte Lysate (LAL)。本研究根據 USP <85>、USP <1085>、USP <1225> 以及 ICH Q2 指南設計了 PSNG rCR 的確效方法，是第一篇直接比較 PSNG rCR 和 LAL 試劑在檢測檸檬酸鈉注射藥物 (Sodium Citrate for Injection) 內毒素時的效能並且確認 PSNG rCR 能有效檢測注射藥物內毒素的研究。

驗證方法為三位分析人員分別在數天內使用 ACC 試劑：PSNG rCR、濁度法試劑 Pyrotell®-T LAL、呈色法試劑 Pyrochrome® LAL 各三個批號，檢測三個批號的檸檬酸鈉注射藥物，每次試驗皆為三重複，再依據各方法學標準線所計算出的內毒素濃度，比較分析三種試劑的確效特性 (Validation Characteristics)：線性 (Linearity)、準確度 (Accuracy)、精密度 (Precision)、範圍 (Range)、最低定量濃度 (Quantitation Limit)、專一性 (Specificity) 和耐變性 (Robustness)，其實驗結果為 PSNG rCR、Pyrotell®-T LAL、Pyrochrome® LAL 皆符合允收標準，也代表 PSNG rCR 和 LAL 具有相同的分析性能 (Analytical performance) (表一、表二)。

Table 3. Assessment of PyroSmart NextGen® Method Suitability Compared to Pyrochrome® According to USP <1085>, USP <85>, USP <1225> and the ICH Q2 Guideline (quantitative test for impurities)

Method Suitability Characteristics	PyroSmart NextGen® Results	Pyrochrome® Results	Acceptance Criteria	
1. Linearity (absolute value, correlation coefficient)	0.05-5.0 EU/mL	0.990	0.994	$ r \geq 0.980$
2. Accuracy (PPC recovery)	Spiked Sample 0.05 EU/mL 0.5 EU/mL 5.0 EU/mL	Min-Max (%) 91-138 122-177 88-134	Min-Max (%) 90-122 127-186 115-152	50-200%
3. Precision 3-1 Repeatability (CV)	Spiked Sample 0.05 EU/mL 0.5 EU/mL 5.0 EU/mL	Min-Max (%) 1-12 7-17 11-19	Min-Max (%) 1-9 7-18 10-23	CV ≤30%
4. Range	0.05-5.0 EU/mL	0.05-5.0 EU/mL	0.05-5.0 EU/mL	Precision, accuracy, and linearity at suitable level
5. Quantitation Limit		At 0.05 EU/mL Accuracy: 91-138% Repeatability: 1-12%	At 0.05 EU/mL Accuracy: 90-122% Repeatability: 1-9%	The lowest concentration of endotoxin that can be quantitatively determined with suitable precision and accuracy
6. Specificity		Lot 1-3 Samples Sample Concentration: <1.11 EU/mL Repeatability: 1-19% PPC Recovery: 88-177%	Lot 1-3 Samples Sample Concentration: <1.11 EU/mL Repeatability: 1-23% PPC Recovery: 90-186%	For a sample matrix that does not contain endotoxin, the endotoxin concentration is determined as undetected with suitable precision and accuracy (PPC recovery)
7. Robustness 7-1 Intermediate Precision (95% CI for CV)	Spiked Sample 0.05 EU/mL 0.5 EU/mL 5.0 EU/mL	Min-Max (%) 9-13 13-18 15-20	Min-Max (%) 8-11 13-18 15-20	CV ≤30%

Note: The pH of the mixture of the reagent and sample solution was determined to be between 6.0 and 8.0.

表一 PSNG rCR 和 Pyrochrome® LAL 具有相同的分析性能

Table 4. Assessment of PyroSmart NextGen® Method Suitability Compared to Pyrotell®-T According to USP <1085>, USP <85>, USP <1225> and the ICH Q2 Guideline (quantitative test for impurities))

Method Suitability Characteristics	PyroSmart NextGen® Results	Pyrotell®-T Results	Acceptance Criteria	
1. Linearity (absolute value, correlation coefficient)	0.003-0.3 EU/mL	0.994	0.996	$ r \geq 0.980$
2. Accuracy (PPC recovery)	Spiked Sample 0.003 EU/mL 0.03 EU/mL 0.3 EU/mL	Min-Max (%) 82-128 107-171 94-129	Min-Max (%) 97-168 158-199 110-133	50-200%
3. Precision 3-1 Repeatability (CV)	Spiked Sample 0.003 EU/mL 0.03 EU/mL 0.3 EU/mL	Min-Max (%) 3-21 2-17 1-11	Min-Max (%) 1-27 1-17 3-14	CV ≤30%
4. Range	0.003-0.3 EU/mL	0.003-0.3 EU/mL	0.003-0.3 EU/mL	Precision, accuracy, and linearity at suitable level
5. Quantitation Limit		At 0.003 EU/mL Accuracy: 82-128% Repeatability: 3-21%	At 0.003 EU/mL Accuracy: 97-168% Repeatability: 1-27%	The lowest concentration of endotoxin that can be quantitatively determined with suitable precision and accuracy
6. Specificity		Lot 1-3 Samples Sample Concentration: <1.11 EU/mL Repeatability: 1-21% PPC Recovery: 82-171%	Lot 1-3 Samples Sample Concentration: <1.11 EU/mL Repeatability: 1-27% PPC Recovery: 97-199%	For a sample matrix that does not contain endotoxin, the endotoxin concentration is determined as undetected with suitable precision and accuracy (PPC recovery)
7. Robustness 7-1 Intermediate Precision (95% CI for CV)	Spiked Sample 0.003 EU/mL 0.03 EU/mL 0.3 EU/mL	Min-Max (%) 13-18 14-19 8-11	Min-Max (%) 13-18 7-10 7-10	CV ≤30%

Note: The pH of the mixture of the reagent and sample solution was determined to be between 6.0 and 8.0.

表二 PSNG rCR 和 Pyrotell®-T LAL 具有相同的分析性能

當原本使用 LAL 的用戶想要導入 PSNG rCR，本研究的比較分析確效方法可以做為應用於自身產品的範例。

原廠完整文章：

[Madeline Kelley, etc., A Demonstration of the Validation Process for Alternative Endotoxin Testing Methods Using PyroSmart NextGen® Recombinant Cascade Reagent, BPB Reports, Vol. 6 No. 2 p.68-75, 2023](#)

PSNG介紹：

[PyroSmart NextGen™-你還在用鰲魚試劑偵測內毒素嗎](#)

岑祥技術部 編撰